# Basic

## SURGICAL SUCTION PUMP

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- FI Käyttöohjeet
- NO Bruksanvisning
- ES Instrucciones de uso
- PT Instruções de utilização

# Table of Contents

<b>Congratulations</b>	<b>4</b>
<b>1 Warnings and safety instructions</b>	<b>5</b>
Safety instructions	8
<b>2 Power supply</b>	<b>9</b>
Plugging in the pump	9
<b>3 Description</b>	<b>10</b>
Introduction	10
Intended use/indications	10
Intended user	10
Intended patient population	10
Important note	10
<b>4 Overview</b>	<b>11</b>
Definition of vacuum	11
Versions and main elements of the pump	11
Back of device	12
Operating elements and indicators	12
Safety Set	12
<b>5 Installation</b>	<b>13</b>
Check initial delivery	13
Remove transport lock	14
Set up mobile version (if available)	14
Assembly of the optional clampholder	14
Set up the Safety Set	14
<b>6 Preparation for use</b>	<b>15</b>
Checks before use	15
Assembly of the basic configuration	15
Assembly of the optional foot switch	16
Assembly of the REUSABLE collection system	16
Assembly of the DISPOSABLE collection system	17
Install equipotential bonding	17
<b>7 Operating instructions</b>	<b>18</b>
Connect Basic to mains power	18
Check maximal vacuum for functional check	18
Changing vacuum level	19
Placing out of operation after use	19
<b>8 Vacuum Assisted Delivery setup</b>	<b>20</b>
<b>9 Changing reusable jar</b>	<b>21</b>

<b>10</b>	<b>Changing disposable liner</b> .....	<b>22</b>
<b>11</b>	<b>Troubleshooting</b> .....	<b>24</b>
	No LED lit .....	24
	Yellow LED indicator lit .....	24
	Motor not running .....	24
	Insufficient vacuum.....	24
<b>12</b>	<b>Replacing defective fuse</b> .....	<b>25</b>
<b>13</b>	<b>Filter test</b> .....	<b>26</b>
<b>14</b>	<b>Cleaning guidelines</b> .....	<b>27</b>
	General notes .....	27
	Medela recommended surface cleaning agents for pump housing .....	27
	Water .....	27
	Cleaning/disinfection machines .....	27
	Disposable products .....	27
<b>15</b>	<b>Warranty and servicing</b> .....	<b>29</b>
	Warranty .....	29
	Servicing/Routine check.....	29
<b>16</b>	<b>Disposal</b> .....	<b>29</b>
<b>17</b>	<b>Accessories overview</b> .....	<b>30</b>
<b>18</b>	<b>Technical specifications</b> .....	<b>31</b>
<b>19</b>	<b>Signs and symbols</b> .....	<b>32</b>
<b>I</b>	<b>Technical documentation (EMC)</b> .....	<b>326</b>
<b>II</b>	<b>List of accessories</b> .....	<b>330</b>

# Congratulations

By choosing the Basic, you have acquired a high-quality suction pump which is adaptable to your needs. As all Medela suction pumps, this pump provides reliable and simple suction. Its simple handling and cleaning as well as the safety features are additional advantages you receive. A comprehensive range of accessories makes the Basic ideally suited to a wide range of medical applications and can be used for continuous operation. Contact us – we will be pleased to advise you.

# 1 Warnings and safety instructions

## WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

## CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

## Safety related tip

Indicating useful information about the safe use of the device.

The Basic is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when the Basic is used in combination with the original Medela accessories (collection system, tubings, filters etc. – see annex “II List of accessories”).

The Basic is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according to clause 7 and 8.9 and can be used in the vicinity of other EMC-tested devices that fulfil the requirements as outlined in the IEC 60601-1-2 standard. HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Basic pump.

**Please read and observe these warning and safety instructions before operation. These Instructions for Use must be kept with the device for later reference.**

Please note that these instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician. Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Basic if it is used in accordance with the Instructions for Use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Subject to change.

## WARNINGS

- For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.
- To avoid risk of electric shock, this equipment must only be connected to a fixed mains socket with protective earth.
- The device must not be used for suctioning explosive, easily flammable or corrosive liquids.
- The connecting tubing supplied with the device must never come into direct contact with the suction area. A sterile suction catheter must always be used (risk of infection).
- Before cleaning the device, pull the plug out of the fixed mains socket.
- No modification of this equipment is allowed.
- Consult the indications for use and consider risk factors and contraindications before using the Basic. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.
- Do not connect this device to a passive drainage tube.
- Not suitable for setting at a low vacuum, as needed for example for thoracic drainage without specialized accessories. Not approved for outdoor use or transport applications.
- The Basic pump may shortly shut down with electrostatic discharge (ESD) events at the DC port of 15kV.

## CAUTIONS

- Incorrect use can cause pain and injury to the patient.
- Do not use sterile accessories when the sterile packaging is damaged.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie can affect the Basic pump and should be kept at least a distance 1ft (30 cm) away from the equipment.
- The rack version requires a minimum distance of 5 cm to the enclosure to prevent overheating of the device.
- The patient should be monitored regularly according to the physicians' instructions and facility guidelines. Objective indications or signs of a possible infection or complication must be met immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge).  
Non-observance can lead to considerable danger of the patient. Monitor the Basic frequently for operating status.
- To prevent the device from overheating, the exhaust at the bottom of the unit must be unobstructed when the unit is operational.

**Safety related tip**

- The Basic Suction Pump is Magnetic Resonance (MR) Unsafe. Do not take the pump into the MR environment.
- For safety tests, the device requires service and repair throughout its service life in accordance with the service manual.
- The protection of the Medela Basic against the effects of the discharge of a cardiac defibrillator is dependent upon the use of appropriate cables.
- Separation from the mains is only assured through the disconnection of the mains plug and the fixed mains socket.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.

## **Safety instructions**

- The Medela Basic is a medical device that requires special safety measures in regard to EMC. It must be installed and put into operation in accordance with the EMC information in annex “I Technical documentation (EMC).”
- In the case of overflow, inform the internal technical service immediately and perform the tasks in the service manual.
- In each of the following cases, the device must not be used and it must be repaired by Medela Customer Service:
  - if the power cord or the plug are damaged
  - if the device is not functioning perfectly
  - if the device is damaged
  - if the device shows clear safety defects.
- Keep the power supply cord away from hot surfaces.
- The mains plug must not come into contact with moisture.
- Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord!
- Never leave the device unattended when it is switched on.
- The pump must stand upright during use.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place the device in water or other liquids.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.

These instructions for use must be kept for later reference.

## 2 Power supply

The Basic is a mains-powered suction pump. Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.

### **Plugging in the pump**

Take the power supply cord and plug into the appliance inlet port at the back of the device. Use the mounting bracket to secure the cord in the inlet port.

# 3 Description

## Introduction

The Basic is a high-quality suction pump. It is powered by a piston/cylinder system and provides maximum suction performance for many suctioning needs. It ideally combines easy handling and cleaning with safety features to ensure optimal operation. You can choose from a comprehensive range of accessories from Medela to configure the pump to many medical applications. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

## Intended use/indications

Intended Use of the Basic suction pump is the creation of a constant vacuum in the range of 0 to -90 kPa. This vacuum can be used for all applications in the hospital, clinic and doctors practice which require a vacuum in this setting such as: general surgery, liposuction, endoscopy, epicardial ablation, nasopharyngeal suction, neurosurgery, OPCAB, vacuum assisted cesarian/delivery and wound drainage.

## Intended user

The Basic should only be operated by properly trained staff. These persons must not be hard of hearing or deaf and must have adequate visual faculty. The training should be refreshed at least once a year.

## Intended patient population

The Basic is intended to be used on patients only exhibiting conditions as described in the indications for use.

## Important note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.


# 4 Overview

## Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in Kilopascal (kPa). In this document, the indication of -10kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079:1999).

## Versions and main elements of the pump


Portable version:


Rack version:

Mobile version:

**Back of device**

**Operating elements and indicators**


**green light** Pump is plugged in

**yellow light** Pump has an error. Refer to chapter 11

**white light** Pump is running

**Safety Set**

# 5 Installation

## 5.1 Check initial delivery

Check the delivery package of the Basic for completeness and general condition.

Basic portable version

**REF** 071.0000

or

Basic rack version

**REF** 071.0001

Mains plug/Allen key

**REF** See service manual

Safety Set

**REF** 077.0711

Silicone Tubing  $\varnothing 7 \times 12$  mm with 2 coupling pieces

**REF** 077.0922

Instructions for use

**REF** 200.5524

## 5.2 Remove transport lock


1.

- 1.1 Remove the red note
- 1.2 Remove 3 screws and store them for later use.


## 5.3 Set up mobile version (if available)

1.

- 1.1 Position top part of trolley on bottom part, making sure the tubing fits as shown.
- 1.2 Connect parts with 4 screws.


2.

- 2.1 Position pump on trolley. Make sure that the front of the pump and the standard rail point forward.
- 2.2 Connect pump with 4 screws.
- 2.3 Attach the cable holder with the two screws to the pump.

## 5.4 Assembly of the optional clampholder (when using the trolley)

1.

- 1.1 Press and hold the blue release knob.
- 1.2 Attach the clampholder to the standard rail by releasing the blue knob.

## 5.5 Set up the Safety Set

1.

- 1.1 Attach the mechanical overflow protection to the lid. Pull gently downwards to make sure it is open/deactivated.
- 1.2 Attach the lid to the jar.
- 1.3 Close the two lid clamps.

2.

- 2.1 Attach the Safety Set to the pump.

# 6 Preparation for use

## WARNINGS

For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.

## CAUTIONS

- The Basic must remain in an upright position during use.
- The rack version requires a minimum distance of 5 cm to the enclosure to prevent from overheating of the device. The back of the enclosure must be open.
- Sterile accessories must be checked to ensure the integrity of the packaging before use.
- Non sterile and reusable accessories must be cleaned, disinfected and/or sterilized according to the Medela cleaning guide (product code 200.2391).

### 6.1 Checks before use

- Check the Basic system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check for completeness and general condition of the Basic delivery package.
- Check all accessories prior to use:
  - suction jars, lids and liners for cracks, brittle and flawed spots. Replace if necessary.
  - tubing for cracks, brittle areas and that connectors are firmly attached. Replace if necessary.
  - as an additional safety test, evacuate the system (including jars) to maximum vacuum before actual use.

### 6.2 Assembly of the basic configuration

1.

1.1 Make sure that the Safety Set is attached to the Basic. See chapter 5.5.


2.

2.1 Attach a filter (if required) to the Safety Set with the arrow pointing in the flow direction.

3. Attach all necessary accessories according to your needs. See annex "II List of accessories".

### 6.3 Assembly of the optional foot switch

- 1.1 Connect the optional foot switch to the pump by plugging in the plug.
- 1.2 Test the correct functioning of the foot switch.

### 6.4a Assembly of the REUSABLE collection system

positions for jars:


#### 1. If you use a DISPOSABLE collection system proceed with step 6.4b

- 2.1 Attach mechanical overflow protection to lid.
- 2.2 Pull gently downwards to make sure it is open/deactivated.

- 3.1 Attach the lid to the jar and lock with the two lid clamps.
- 3.2 Attach the jar to the carrying handle (portable version only) or to the optional clampholder on the rail (see step 5.4).

- 4.1 Connect tubing from the Safety Set to the lid of the jar (vacuum port).

- 5.1 Connect patient tubing to the lid of the jar (patient port).

## 6.4b Assembly of the DISPOSABLE collection system

1. If you use a REUSABLE collection system proceed with step 6.4a


2.1 Choose the size necessary for the application.  
Yellow colour code: 1.5l,  
Orange colour code: 2.5l.


2.2 Attach the jar to the carrying handle (portable version only) or to the optional clampholder on the trolley (see step 5.4).

3.1 Prepare and insert corresponding liner size according to instructions delivered with liners. The disposable liner has an integrated overflow protection filter. No action is necessary to activate it.

4.1 Connect tubing from the jar (A) to the Safety Set (B).


5.1 Connect patient tubing to the lid of the liner using either the coloured angle piece or the patient port directly (depending on tubing size).

## 6.5 Install equipotential bonding

to bus bar

and

1.1 Connect the equipotential conductor at the back of the Basic with the bus bar as advised by your technician and according to inhouse directives. The respective cable is not supplied from Medela.

The equipotential conductor and connecting cables are used to equalise the potential between equipment placed in the hospital environment.

# 7 Operating instructions


## CAUTIONS

The Basic is to be set up in such a way, that a separation from the mains supply can be easily managed.

### 7.1 Connect Basic to mains power

1. Check the pump before use following the instruction in chapter 6.1.

2.1 Connect the power cord to the appliance inlet at the back of the Basic.  
2.2 Plug in the mains plug of the power cord to a fixed mains socket.

3.1 An internal self-test is performed. When the green LED lights up, the device is ready for use.

### 7.2 Check maximal vacuum for functional check

1.1 Switch on the Basic.

2.1 Turn the vacuum regulator to the right to set maximum vacuum

3.1 Seal the end of the patient tubing with your thumb.  
3.2 Compare the maximum vacuum according to the specification (on left). See chapter 11 if the maximum vacuum is not reached.

#### Specifications:

Altitude  
above  
sea level:

Max. Vacuum:

+ 2000 m	- 70 kPa - 525 mmHg
+ 1000 m	- 79 kPa - 592 mmHg
+ 500 m	- 85 kPa - 638 mmHg
0 m	- 90 kPa - 675 mmHg

(Tolerance: +/- 15%)

**CAUTIONS**

When the Basic is used for wound drainage, the negative pressure should be set according to instruction of the specialist and not cause any wound damage.

**7.3 Changing vacuum level**


- 1.1 Clamp patient tubing
- 1.2 Turn vacuum regulator to select the correct vacuum according to the particular application.
- 1.3 Check vacuum gauge for setting.

**7.4 Placing out of operation after use**

- 1.1 Touch on/off button to switch off the Basic.

- 2.1 Disconnect the mains plug from the fixed mains socket.

- 3. Clean and disinfect the Basic. See chapter 14.

# 8 Vacuum Assisted Delivery setup


## CAUTIONS

The Basic is to be set up in such a way, that a separation from the mains supply can be easily managed.

1. Set up pump and accessories according to chapter 6.


2.

2.1 Connect the foot vacuum regulator: the silver adapter must be fully seated in the Safety Set of the suction pump.

2.2 Attach tubing to top of metal adapter.

3a

3.1 Attach tubing from suction cup to patient connection on the lid of the liner.

or

3b

3.1 Attach tubing from suction cup to patient connector on suction jar of the reusable collection system.

## Specifications:

Altitude

above

sea level:

Max. Vacuum:

+ 2000m	- 70 kPa - 525 mmHg
+ 1000m	- 79 kPa - 592 mmHg
+ 500m	- 85 kPa - 638 mmHg
0m	- 90 kPa - 675 mmHg

(Tolerance: +/- 15%)

4.

4.1 Switch on pump, turn to max. vacuum, clamp tubing from suction cup and fully depress the vacuum regulator (forward and down, using ball of foot).

4.2 Compare maximum vacuum according to specification (on left).


5.

5.1 If OK, release vacuum by returning foot vacuum regulator to resting «zero pressure» state (rear and down using heel of foot).

5.2 The pump is now ready for use.

# 9 Changing reusable jar

- 1.1 Switch pump off by touching the on/off button.
- 1.2 Remove patient and vacuum tubings from the lid.
- 1.3 Remove full jar from the Basic.

- 2.1 Prepare new jar and unpack it.
- 2.2 Secure new jar in the clampholder or directly to the Basic carrying handle.
- 2.3 Reconnect vacuum tubing and new patient tubing firmly.

3. Empty jar and dispose of patient tubing according in accordance with local guidelines and inhouse guidelines.

- 4.1 Switch pump on by touching the on/off button. Vacuum is built up.
- 4.2 Set vacuum according to the particular application. See chapter 7.3.

# 10 Changing disposable liner


## WARNINGS

Disposable liners are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics.

1.

1.1 Remove patient tubing with coloured angle piece from lid.

1.2 Plug patient port.

2.

2.1 Switch pump off by touching the on/off button.

3.

3.1 Reduce vacuum to minimum by turning vacuum regulator to the left.

4.


4.2 Remove full liner from the jar.

5. Dispose of the liner and patient tubing in accordance with local guidelines and inhouse guidelines.

6.

- 6.1 Prepare new liner and unpack it.
- 6.2 Insert new liner in the jar and press down lid.

7.

- 7.1 Connect tubing from the jar to the Safety Set.
- 7.2 Connect new patient tubing firmly.

8.

- 8.1 Switch pump on by touching the on/off button. Vacuum is built up.
- 8.2 Set vacuum according to the particular application. See chapter 7.3.

# 11 Troubleshooting

## **No LED lit**

The Basic is not connected to the mains or the fuse needs replacement.

## **Yellow LED indicator lit**

### **Minor case**

yellow LED indicator lit but the pump can be switched on and off:

- contact the internal technical department or your authorised service center at next possible occasion.

### **Major case**

yellow LED indicator lit and pump cannot be switched on and off:

- contact the internal technical department or your authorised service center for repairs / maintenance.

## **Motor not running**

Check to ensure:

- the Basic is switched on. The standby LED must be illuminated.
- the mains plug is inserted correctly into the fixed mains socket and into the appliance inlet.
- the fuse on the back of the Basic is not defective. For replacing the defective fuse see chapter 12.

If the fault cannot be rectified, please contact the internal technical department.

## **Insufficient vacuum**

Check if:

- the vacuum regulator is set correctly.
- the tubing is not defective or broken. If necessary, replace.
- all plug-in connections are tight.
- the overflow protection is deactivated / open. If the overflow protection is activated, deactivate it as shown under 6.4a/2.2.
- the suction jar / lid have no cracks, brittle areas, discolouration. If necessary, replace.
- the disposable system has no cracks, brittle areas, discolouration. If necessary, replace.
- the filter is not clogged. To test if the filter is clogged see chapter 13.

If the fault cannot be rectified, please contact the internal technical department.



# 12 Replacing defective fuse

## WARNINGS

Before replacing the fuse, pull the mains plug from Basic out of the fixed mains socket.

1. 
  - 1.1 Disconnect the Basic from the mains.
  - 1.2 If attached, remove the cable holder (2 screws on bottom).
  
2. 
  - 2.1 Open (pull out) the fuse holder on the back of the Basic.
  
3. Replace the defective fuse(s). Make sure that the technical specifications are correct.  
**T 1.6AH, 250VAC, 5x20mm**
  
4. 
  - 4.1 Close (push back) the fuse holder.
  
5. 
  - 5.1 Reconnect the Basic to the mains again.
  - 5.2 Switch on the Basic again.

# 13 Filter test

1.

1.1 Remove tubing from filter.

2.

2.1 Select maximum vacuum.

3.

3.1 Switch on Basic by pressing the on/off button.

4.

4.1 Read the vacuum.  
4.2 Replace the filter if the vacuum exceeds -20kPa.

# 14 Cleaning guidelines

## WARNINGS

After each use, the parts that have been in contact with the aspirated secretions are to be cleaned, disinfected, sterilised or disposed of **according to the table on the next page**.

**Before cleaning the device, pull the mains plug out of the fixed mains socket.**

## Safety related tip

Avoid contact of fluids with the ends of the mains plug or appliance inlet port.

## General notes

- These are general recommendations only that may be adjusted individually, based on the hospital's specific directives and cleaning practices and policies.
- Also refer to detailed information according to Medela cleaning instructions (product code 200.2391).
- Wear protective gloves for cleaning/disinfection.
- Dispose of fluids such as blood and secretions and the parts contaminated with them in accordance with internal hospital guidelines.

## Medela recommended surface cleaning agents for pump housing

- **Perform** [concentrate], Schülke & Mayr, [www.schuelkemayr.com](http://www.schuelkemayr.com)

## Water

Use only the purest quality of water for cleaning. Water hardness is a serious consideration since deposits left on medical products may not be properly decontaminated. Use deionised water in order to reduce this problem. The final rinse water should be bacterial free and contains no endotoxins.

## Cleaning/disinfection machines

Can be used to desinfect parts from the table on the next page. A hot water rinse (maximum temperature 100 °C) may provide a medium-to high level of disinfection. Every section of the constituent parts must be accessible in order to ensure efficient cleaning. We recommend using a cleaning/disinfecting machine that has been approved by the Robert Koch-Institute and complies with ISO 15883. Recommended temperature for noncritical medical devices (i.e. those that only come into contact with uninjured skin) is 90 °C for 1 minute. The time is increased to 5 minutes for all medical devices that are considered to be critical.

## Disposable products

These are single use products not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

PSU* jars, lids	PC** jars	Silicone tubing	Reusable suction cups	Pump housing Mainscable	Plastic accessories	Legend: * Polysulfone ** Polycarbonate
x	x	x	x	x	x	<b>1. Disassembly</b> Separate all individual parts before cleaning, disinfecting and sterilising.
x	x	x	x			<b>2. Cleaning</b> Clean components in hot water (60–70 °C) containing a detergent with a pH range between 6.0 and 8.0 only, in order to avoid damaging the instruments and containment devices. Enzymatic detergents help with the removal of organic matter, such as blood. Detergents should be used according to their manufacturer's recommended concentration levels. Some alkaline detergents have been formulated to be safe for reprocessing medical devices. The manufacturers of these detergents should provide information about specific materials that may be damaged by their detergents. Soak all parts thoroughly with warm, soapy water or in enzymatic detergent for 1–5 minutes. <b>1.</b> Remove visible dirt with a cleaning tool – general purpose cleaning brushes, such as pipe cleaners or non-abrasive lint cloths. Brushes and pipe cleaners should fit snugly but still be able to be moved around easily within the area to be cleaned. Rinse thoroughly in clear water. Allow to dry. <b>2.</b> Check the instruments for visible dirt and repeat these steps if necessary.
				x	x	<b>3. Cleaning</b> Wipe with detergent mentioned above. Use clean cloth to dry.
x	x	x	x		x	<b>4. Disinfect</b> Soak at room temperature for 30 minutes in a disinfection/cleaning solution (A). After the reaction time, use water to rinse the solution residue from the individual parts and rinse the tubing. Rinse for at least 2 minutes with the cleaning solution (B).  (A) fluid aldehyde-free disinfection solution with cleaning effect for use as a soaking bath, bactericide, fungicide, limited virucide, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone, slightly alkaline. e.g. 2% neodisher® Septo MEDsolution  (B) Fluid, pH-neutral, enzymatic cleaning solution for the treatment of instruments by machine or manually with very good material compatibility for stainless steel, non-ferrous material and plastics including silicone.  or use cleaning/disinfecting machine.
x		x	x			<b>5. Sterilise</b> Remove any chemical disinfectant residue before autoclaving. Do not stack during autoclaving. In the autoclave at 134 °C for 10 minutes. The sterilisation device must comply to ISO 17665:2003. The following procedures are acceptable: In the autoclave with saturated steam at 132 °C, with triple pre-evacuation, sterilisation time of 6 minutes, or at 135–137 °C in an porous load vacuum autoclave for 3–3.5 minutes
x	x	x	x		x	<b>6. Storing components</b> After disinfection/sterilising, store the components in sterile foil until required for use.

## 15 Warranty and servicing

### Warranty

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 5 years from the date of delivery. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our pumps. In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on the pump by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

### Servicing/Routine check

Routine checks and service work are only to be carried out by positions authorised by Medela. Medela recommends to carry out the routine check 1x per year (see service manual in English). The Medela service manual, is available upon request.

## 16 Disposal

The Basic is made of metals and plastics and should be disposed of in accordance with the European directives 2011/65/EU and 2012/19/EU. Additional, local guidelines must also be observed. Please take care that you dispose of the Basic and its accessories in accordance with your local disposal guidelines.

### User information for the disposal of electrical and electronic equipment

This symbol means that the electrical and electronic equipment must not be disposed together with unsorted municipal waste. Correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

# 17 Accessories overview

## WARNINGS

Basic was verified in combination with the accessories listed in annex “II List of accessories”. For correct and safe operation, use Basic with these accessories only. Further information is supplied on the instruction sheet of the individual accessory.

## Safety related tip

If the pump is used together with third party patient interfacing devices (e.g. cannulas, catheters), they must:

- have CE mark (within the European Union) and local registration (outside the European Union)
- be able to be attached to Medela accessories safely without impacting the performance of the pump.

    Tubing connection on reusable lids of jars: Ø6–10mm, Ø10–14 mm

    Tubing connection on disposable liners: Ø6.5–11 mm

**Note:** when combining Medela parts and a new patient interfacing device, you take over the responsibility of the entire system and should test the combination to ensure the vacuum levels are properly maintained.

# 18 Technical specifications

high vacuum  
- 90 kPa / - 675 mmHg  
Tolerance: +/- 15 %

Measured at 0m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

high flow  
30l/min.  
+/- 5l/min.

9.3 kg  
20.5 lbs  
Rack version


100-240V, 50/60Hz  
50 W


ISO 13485  
CE (93/42/EEC), Ila

HxWxD (rack version)  
210x305x375 mm  
8.2x 12.0x 14.8 inches

Transport/Storage Conditions

Operating Conditions

IP21

# 19 Signs and symbols

**0123** This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the class of the pump.

This symbol indicates the date of manufacture.

This symbol indicates protective earth (ground).

This symbol indicates that the device should not be used after the date shown.

**IP21**

This Symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

**REF**

This symbol indicates manufacturer's catalogue number.

This symbol indicates the legal specifications of the pump.

This symbol indicates a type CF applied part.

**SN**

This symbol indicates manufacturer's serial number.

This symbol indicates to follow instructions for use.

This symbol indicates MR UNSAFE.

**LOT**

This symbol indicates manufacturer's batch code.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates that interferences may occur in the vicinity of equipment marked with this symbol.

**STERILEEO**

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates the equipotential conductor connection.

This symbol indicates safety related tip.

This symbol indicates the connection port for the foot switch.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates earth (ground).

This symbol indicates the manufacturer.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates to handle the fragile device with care.

This symbol indicates the fuse

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates do not use the device if package is damaged.

This symbol indicates the maximum vacuum level of the pump.

This symbol indicates the flow levels of the pump.

This symbol indicates the number of items n that the content is sufficient for.

This symbol indicates the electrical specifications of the pump.

This symbol indicates to not dispose the device together with unsorted municipal waste (for EU only).

This symbol indicates the weight of the pump.

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates the dimensions (h x w x d) of the pump.

This symbol indicates a carton package.

This Symbol indicates a Prescription Device. U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).

This symbol indicates to keep the device away from sunlight.

This symbol indicates number of items.

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# I Technical documentation

## WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Basic pump.

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Basic pump.

## EMC (Electromagnetic Compatibility)

Basic has been EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according to clause 7 and 8.9 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements of the relevant IEC 60601-1-2 standard. Basic is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the following EMC information.

## Essential Performance

Medela has determined that the Basic suction pump does not have any functions, the failure of which to work properly, could impact the safety of the system. As such, there is no essential performance as defined in the IEC 60601-1-2 standard. Nevertheless, the Basic was tested in accordance with IEC 60601-1-2 ensuring that the pump could maintain the correct vacuum during testing.

Electromagnetic Compatibility (EMC, IEC 60601-1-2, 3rd Edition and IEC60601-1-2 4th Edition, Table 1)

## Electromagnetic Emissions

Basic is intended for use in the electromagnetic environment specified below. The customer or the user of Basic should assure that it is used in such environment.

Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Basic uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Basic is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 60000-3-3	Complies	

## WARNINGS

Basic should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Basic should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Compatibility (EMC, IEC 60601-1-2, 3rd Edition, Table 2)

## Electromagnetic Immunity

Basic is intended for use in the electromagnetic environment specified below. The customer or the user of Basic should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70% $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95 % dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70% $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Basic requires continued operation during power mains interruptions, it is recommended that Basic be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	30A/m	Power frequency magnetic fields should be at levels of a typical location in a typical commercial or hospital environment.

**NOTE**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Electromagnetic Immunity

Basic is intended for use in the electromagnetic environment specified below. The customer or the user of Basic should assure that Basic is used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Basic, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{P}$ $d = 0.35 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 0.7 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Basic used exceeds the applicable RF compliance level above, Basic should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Basic.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

### Recommended separation distance between portable and mobile RF communications equipment and Basic

Basic is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Basic can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Basic as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
W			
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.11	1.11	2.21
100	3.5	3.5	7


For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# II List of accessories

### Patient tubing (from collection system to patient applied part)

<p><b>Disposable, 7x10 mm, sterile</b></p>  077.0170: 150 cm	<p><b>Silicone, 7x12 mm</b></p>  077.0053: 32 cm 077.0055: 120 cm 077.0056: 150 cm 077.0051: 200 cm 077.0164: 300 cm 077.0052: 2500 cm
<p><b>Disposable, 6x8 mm, sterile</b></p>  077.0184: 150 cm	<p> 077.0961: 180 cm, silicone (green)</p>
<p><b>Disposable, 6x8 mm</b></p>  077.0951: 180 cm, sterile  077.0952: 180 cm  077.0191: 300 cm, sterile  077.0193: 300 cm OEM, sterile	<p><b>Disposable, 12x18 mm</b></p>  077.0963: 200 cm
	<p><b>Silicone, 5x10 mm</b></p>  077.0970: 2500 cm
	<p> 077.1017</p> <p> 077.1020</p>

\*

#### Liners

077.0083/84 077.0086/87  
 077.0194 OEM/ 101033145/  
 101033142/ 101033143  
 101032951/  
 101035340 OEM

#### Reusable lids

077.0450 077.0440 077.0420/30

#### Reusable sets

077.0701: 1l  
 077.0702: 2l  
 077.0703: 3l  
 077.0705: 5l  
 077.0715: 1l  
 077.0716: 2l

#### PC jars (max. 4x2.5l)

077.0082 077.0085

#### Reusable jars (max. 4x5l)

077.0110: 1l 077.0125: 0.25l  
 077.0120: 2l 077.0155: 0.5l  
 077.0130: 3l  
 077.0150: 5l



077.0704: 2x3l  
 077.0706: 2x5l  
 077.0707: 2x2l

### Vacuum tubing (from suction pump to collection system)

**Vacuum tubing, silicone, 7x12 mm**




 077.0054: 50cm	 077.0921: 30cm	 077.1019
 077.0901: 100cm	 077.0922: 60cm	
 077.0902: 200cm	 077.0189: 300cm	
 077.0905: 500cm		
 077.0900: 2500cm	 077.0154: 25cm	 077.1018
 077.0931: 25cm	 077.0911: 33cm	
 077.0185: 60cm	 077.0912: 60cm	 077.1015/16
	 077.0913: 100cm	

Filters (optional)						
			077.0572	077.0575	101035264	101035286
077.0581	077.0092	077.0571	077.0577	077.0576	101035265	101035287
			077.0573	077.0582	101035284	101035288
			077.0574	101035263	101035285	101035289

**Tubing for high volume setups**

	077.0563	Change-over valve		077.0158:	34 cm (100 pc) serial connection DCS, disposable: 8x12 mm
	077.0542:	50 cm serial connection RCS, silicone: 7x12 mm		077.0160:	34 cm serial connection DCS, disposable: 8x12 mm
				077.0095:	25 cm, serial connection, vacuum tubing, DCS

- Patient tubing / tubing connectors / Patientenschläuche / Verbindungsstücke / Tuyaux patients / tuyaux de connexion à l'aspirateur / Tubi paziente / raccordi per tubi / Patiëntslangen / slangaansluitstukken / Kopplingar för patientslangar / slangar / Potilasletkusto / letkuston liittimet / Pasientslanger / slangekoblinger / Tubos de paciente / conexiones de tubo / Tubagem do doente / conectores da tubagem
- Liners / Beutel / Poches / Sacche / Opvangzakken / Påsar / Tiivisteet / Poser / Bolsas / Mangas
- Reusable lids / Mehrweg-Deckel / Couvercles réutilisables / Coperchi riutilizzabili / Herbruikbare deksels / Flergångslock / Kestokäyttöiset kannet / Flergangs lokk / Tapas reutilizables / Tampas reutilizáveis
- Reusable sets / Mehrweg-Sets / Sets réutilisables / Kit riutilizzabili / Herbruikbare sets / Flergångssatser / Kestokäyttöiset setit / Flergangs sett / Kits reutilizables / Kits reutilizáveis
- PC jars / PC-Behälter / Bocaux en PC / Contenitori in PC / PC-opvangpotten / PC-behållare / PC-imupurkit / PC-beholdere / Recipientes de PC / Recipientes PC
- Reusable jars / Mehrweg-Behälter / Bocaux réutilisables / Contenitori riutilizzabili / Herbruikbare opvangpotten / Flergångsbehållare / Kestokäyttöiset imukupit / Flergangs beholdere / Recipientes reutilizables / Recipientes reutilizáveis
- Vacuum Tubing / Vakuumschlauch / Tubulure à vide / Tubi per il vuoto / Vacuümslang / Vakuümslang / Tyhjiöletku / Vakuümslanger / Tubos de vacío / Tubos de vacío
- Filters / Filter / Filtres / Filtri / Filters / Filter / Suodattimet / Filtre / Filtros / Filtros
- Foot controls / Fußregler / Pédales / Comandi a pedale / Voetregelaars / Fotkontroller / Jalkasäätimet / Fotkontroller / Interruptores de pie / Controlos de pedal
- Clampholders / Schienenklemmen / Accessoires, si le rail standard est accessible / Portastaffa / Accessoires indien uitgevoerd met standaardrail / Klover / Kiskokiinnittimet / Klemmeholdere / Soportes de pinza / Soportes dos ganchos
- Applied parts Vacuum Assisted Delivery / Anwendungsteile für Vakuumextraktion / Accessoires destinés à l'accouchement par ventouse / Parti applicate per il parto assistito con vuoto / Toegepaste onderdelen voor vacuümbevallingen / Tillbehör för sugklockeförlossning / Käytettävät osat, imukupisynnytyks / Deler for bruk ved vakuumassistert forløsning / Componentes aplicados para el parto asistido con ventosa / Peças aplicadas em parto assistido por vacío

- Specimen cups/Gewebe-Auffangbehälter/Collecteurs de tissus/Coppette per campioni/  
Monsterpotjes/Vävnadsfälla/Näytesäiliöt/Prøvetakingsbegere/Recipiente de muestr-  
ras/Copos para amostra
  
- Drainage valve/Drainageventil/Valve de drainage/Valvola di drenaggio/  
Drainageventiel/Dräneringsventil/Tyhjennysventtiili/Dreneringsventil/Válvula de drenaje/  
Válvula de drenagem
  
- Additional accessories/Weiteres Zubehör/Accessoires supplémentaires/Accessori  
supplementari/Overige accessoires/Ytterligere tillbehør/Muut lisävarusteet/Ytterligere  
tilbehør/Accesorios adicionales/Acessórios adicionais
  
- Tubing for high volume setups/Schlauch für hochvolumige Aufbauten/Tubulure pour  
installations à grand volume/Tubi per configurazioni a volume elevato/Slang voor  
hoogvolume-opstelling/Slang för stora volymer/Suuren tilavuuden asennuksiin  
tarkoitettu letku/Slanger for oppsett med høyt volum/Tubos para configuraciones con  
grandes volúmenes/Tubos para configurações de volume elevado
  
- serial connection, vacuum tubing, Disposable Collection System (DCS)/Serienschaltung,  
Vakuumschlauch, Einweg-Auffangsystem (DCS)/raccord en série, tubulure à vide, système  
de collecte à usage unique (DCS)/collegamento in serie, tubi per il vuoto, sistema di  
raccolta monouso (DCS)/seriële aansluiting, vacuümslang, disposable opvangsysteem  
(DCS)/serieanslutning, vakuümslang, engångssystem för uppsamling av vätska (DCS)/  
sarjaliitântä, tyhjiöletku, kertakäyttöinen keräysjärjestelmä (DCS)/seriekobling,  
vakuümslanger, Oppsamlingssystem for engangsbruk (DCS)/conexión en serie, tubos de  
vacío, sistema de recogida desechable (DCS)/ligação em série, tubos de vácuo, Sistema  
de Recolha Descartável (SRD)
  
- serial connection, Disposable Collection System (DCS), disposable tubing/  
Serienschaltung, Einweg-Auffangsystem (DCS), Einweg-Schlauch/raccord en série,  
système de collecte à usage unique (DCS), tubulure à usage unique/collegamento in serie,  
sistema di raccolta monouso (DCS), tubi monouso/seriële aansluiting, disposable  
opvangsysteem (DCS), disposable slang/serieanslutning, engångssystem för uppsamling  
av vätska (DCS), engångsslang/sarjaliitântä, kertakäyttöinen keräysjärjestelmä (DCS),  
kertakäyttöinen letku/seriekobling, Oppsamlingssystem for engangsbruk (DCS),  
engångsslanger/conexión en serie, sistema de recogida desechable (DCS), tubos  
desechables/ligação em série, Sistema de Recolha Descartável (SRD), tubos descartáveis

– serial connection, Reusable Collection System (RCS), silicone/Serienschaltung, Mehrweg-Auffangsystem (RCS), Silikon/raccord en série, système de collecte réutilisable (RCS), silicone/collegamento in serie, sistema di raccolta riutilizzabile (RCS), silicone/seriële aansluiting, herbruikbaar opvangsysteem (RCS), silicone/serieanslutning, flergångssystem för uppsamling av vätska (RCS), silikon/sarjaliitämä, uudelleenkäytettävä keräysjärjestelmä (RCS), silikon/seriekobling, Oppsamlingsssystem for gjenbruk (RCS), silikon/conexión en serie, sistema de recogida reutilizable (RCS), silicona/ligação em série, Sistema de Recolha Reutilizável (SRR), silicone

Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
www.medelahealthcare.com

#### International Sales

Medela AG  
Lättichstrasse 4b  
6341 Baar  
Switzerland  
Phone +41 41 562 51 51  
Fax +41 41 562 51 00  
customer@medela.ch  
www.medelahealthcare.com

#### Australia

Medela Australia Pty Ltd,  
Medical Technology  
3 Arco Lane  
Heatherston, Vic 3202  
Australia  
Phone +61 3 9552 8600  
Fax +61 3 9552 8699  
contact@medela.com.au  
www.medelahealthcare.com.au

#### Canada

Medela Canada Inc.  
4160 Sladeview Crescent Unit # 8  
Mississauga, Ontario  
L5L 0A1  
Canada  
Phone +1 800 435 8316  
Fax +1 800 995 78 67  
info@medela.ca  
www.medela.ca

#### China

Medela (Beijing)  
Medical Technology Co., Ltd.  
Unit M, 4F  
No 1068 Wuzhong Road,  
Minhang District  
Shanghai 201103  
P.R. of China  
Phone +86 21 6236 8102  
Fax +86 21 6236 8103  
info@medela.cn  
www.medela.cn

#### France

Medela France SAS  
14, rue de la Butte Cordière  
91154 Etampes cedex  
France  
Phone +33 1 69 16 10 30  
Fax +33 1 69 16 10 32  
info@medela.fr  
www.medela.fr

#### Germany & Austria

Medela Medizintechnik  
GmbH & Co. Handels KG  
Postfach 1148  
85378 Eching  
Germany  
Phone +49 89 31 97 59-0  
Fax +49 89 31 97 59 99  
info@medela.de  
www.medela.de

#### India

Medela India Private Limited  
602-603 | 6th Floor,  
JMD Megapolis | Sohna Road,  
Sector - 48 | Gurgaon 122 001  
India  
Phone +91 124 4832069  
Fax +91 124 4832066  
info@medela.in  
www.medela.in

#### Italy

Medela Italia Srl a socio unico  
Via Turrini, 13/15 – Loc. Bargellino  
40012 Calderara di Reno (BO)  
Italy  
Phone +39 051 72 76 88  
Fax +39 051 72 76 89  
info@medela.it  
www.medela.it

#### Japan

Medela K.K.  
G-PLACE Sasazuka 13F  
2-26-2, Sasazuka, Shibuya-ku  
Tokyo 151-0073  
Japan  
Phone +81 3 3373 3454  
Fax +81 3 3373 3457  
info@medela.jp  
www.medela.jp

#### Netherlands & Belgium

Medela Benelux BV  
Uilenwaard 31  
5236 WB's-Hertogenbosch  
The Netherlands  
Phone +31 73 690 40 40  
Fax +31 73 690 40 44  
info@medela.nl  
info@medela.be  
www.medela.nl  
www.medela.be

#### Poland

Medela Polska Sp. z o.o.  
Ul. Lewinowska 8  
03-684 Warszawa  
Poland  
Phone +48 22 865 12 50  
Fax +48 22 864 38 65  
info@medela.pl  
www.medela.pl

#### Russia

OOO Medela  
Vavilova str. 97, entrance 3  
Office 16-17  
117335 Moscow  
Russia  
Phone +7 495 980 6194  
Fax +7 495 980 6194  
info@medela-russia.ru  
www.medela-russia.ru

#### Spain & Portugal

Productos Medicinales Medela, S.L.  
C/ Llacuna, 22 4a Planta  
08005, Barcelona  
Spain  
Phone +34 93 320 59 69  
Fax +34 93 320 55 31  
info@medela.es  
info@medela.pt  
www.medela.es  
www.medela.pt

#### Sweden

Medela Medical AB  
Box 7266  
187 14 Täby  
Sweden  
Phone +46 8 588 03 200  
Fax +46 8 588 03 299  
info@medela.se  
www.medela.se

#### Switzerland

Medela AG  
Lättichstrasse 7  
6341 Baar  
Switzerland  
Phone +41 848 633 352  
Fax +41 41 562 14 42  
contact@medela.ch  
www.medela.ch

#### United Kingdom

Medela UK Ltd.  
Huntsman Drive  
Northbank Industrial Park  
Irlam, Manchester M44 5EG  
UK  
Phone +44 161 776 0400  
Fax +44 161 776 0444  
info@medela.co.uk  
www.medelahealthcare.co.uk

#### USA

Medela LLC  
1101 Corporate Drive  
McHenry, IL 60050  
USA  
Phone +1 877 735 1626  
Fax +1 815 307 8942  
info-healthcare@medela.com  
www.medelahealthcare.us